

# Biacore Tutorial

---

---

## I. Introduction.

The majority of this tutorial is derived from the Biacore manual, *Getting Started BIACORE® 3000*. This tutorial is designed to acquaint the user with the basic operating principles and procedures of the Biacore 3000. In general, most experiments will follow this same format although there will be differences in the coupling step depending on which sensor chip and/or immobilization chemistry is used.

For Further Reading:

1. *Getting Started BIACORE® 3000*
2. *BIACORE 3000® Instrument Handbook*
3. *BIAtchnology Handbook*

## II. The Basics:

Biacore experiments are setup in the following manner generally:

- a. Immobilization: One of the interactants is immobilized to the surface of the sensor chip using coupling chemistry. The immobilized interactant is thus referred to as the **ligand**.
- b. Binding: The interaction partner is injected across the ligand-immobilized surface resulting in binding. The interaction partner in solution is referred to as the **analyte**.
- c. Regeneration: The bound analyte is removed from the surface of the sensor chip through the injection of a regeneration solution.
- d. Controls: Controls should be set up during the experiment to test for non-specific binding and to verify that the observed binding is from the interaction being studied.

The steps covered in this tutorial are specific to the interaction of myoglobin to immobilized anti-myoglobin antibody. While this procedure is designed for this antibody-antigen interaction, the protocol can be easily adapted for another interaction. The following steps will be covered in this tutorial:


- a. Immobilization of Anti-Myoglobin Monoclonal Antibody ( $\alpha$ -Myo MAB) via amine coupling. Optimal immobilization conditions will be sought through a series of test injections.
- b. Myoglobin will be injected and the optimal regeneration conditions for the  $\alpha$ -Myo MAB will be deduced.
- c. Non-specific binding of myoglobin will be tested through injection over a blank reference surface.

## III. Pre-Experiment Notes.

- a. Buffer: Buffer should always be degassed and filtered (0.22um filter) before diluting the samples. The purpose of degassing is to prevent the injection of air into the system which can cause signal spikes that obscure the data. Filtering prevents the injection of particulates into the system which may non-specifically bind and/or clog the microfluidics system.
  - b. Be sure to have all the reagents in the proper rack position before starting a method or using the manual commands. Incorrect rack positions can cause the injection needle to hit the rack and may cause the needle to bend.
  - c. Instrument damage is a serious matter and if it occurs, the user will be billed for the troubleshooting, repairs and replacement parts necessary to get the instrument running again.
  - d. It is **essential** to always leave the instrument with buffer or water running at a low flow rate post-experiment. Letting the instrument run dry can cause serious damage.
  - e. The Biacore supplies (vials, caps, etc.) in the core are part of the per-hour cost of the instrument and are free for use. Sensor chips are available in the facility but are available for purchase at the facility's cost.
  - f. The nomenclature used in Biacore to denote rack positions is as follows:
- 
-

## Biacore Tutorial

2



**Figure 1.** Rack Positions demonstrated. R1 always corresponds to the leftmost reagent rack in the instrument while R2 corresponds to the rightmost

### IV. Materials

1. CM5 Sensor Chip
2. N-(3-Dimethylaminopropyl)-N'-ethylcarbodiimide (EDC)
3. Ethanolamine
4. N-Hydroxysuccinimide (NHS)
5. P20, P200 and P1000 Pipettors and tips
6. Biacore Vials
7. Eppendorf Vials (0.6mL)
8. BiaDesorb Solutions 1 and 2
9. The Following Solutions (All Found in the Getting Started Reagent Kit):

HBS-EP buffer, 200 ml.	10 mM HEPES, pH 7.4, 150 mM NaCl, 3 mM EDTA, 0.005% (v/v) Surfactant P20.
Ligand, 80 $\mu$ l	$\alpha$ -myoglobin mab (monoclonal antibody), 1 mg/ml. The ligand should be used immediately after dilution.
Analyte, 100 $\mu$ l.	Sheep myoglobin, 75 $\mu$ g/ml. The analyte can be stored for one week after dilution.
Coupling buffer 1, 1.5 ml.	5 mM maleate, pH 7.0.
Coupling buffer 2, 1.5 ml.	5 mM maleate, pH 6.0.
Coupling buffer 3, 1.5 ml.	5 mM maleate, pH 4.8.
Regeneration solution 1, 1.5 ml.	10 mM acetate, pH 4.0.
Regeneration solution 2, 1.5 ml.	10 mM glycine, pH 3.0.
Regeneration solution 3, 1.5 ml.	10 mM glycine, pH 2.0.
Green sucrose solution, 1.5 ml.	15% (w/v) sucrose in HBS-EP buffer with color agent.
BIA normalizing solution, 70%, 1.5 ml.	70% (w/w) glycerol.
Empty vials for sample dilution, 18 each.	Polypropylene, 7 mm diameter.
Empty vials for normalization, 6 each.	Glass, 9 mm diameter.

**Fig 2.** Reagent Solutions Necessary for the Tutorial

## Biacore Tutorial

---

---

### V. Preparing the System with a New Sensor Chip.

1. In the Biacore Software, go to the **Tools Menu** and select **Working Tools**. Under working tools select **Desorb** and follow the on-screen instructions. Desorb should be run prior to any to clean the system of adsorbed material. This procedure takes approximately 15mins and should always be run with the Sensor Maintenance Chip docked to the system.
2. Undock the Sensor Maintenance Chip from the system by going to the **Command Menu** and selecting **Undock**. Check **Empty Flow Cells** and hit **Undock**. Wait for the sensor chip to undock, at which point the Sensor Chip light on the front of the instrument will begin to blink.
3. Pull down the sensor chip cover on the front of the instrument, slide the cover back to release the chip and remove the Sensor Maintenance Chip from the instrument. Place the CM5 Sensor Chip in the instrument, lock the sensor chip into place and close the cover on the port.
4. Select **Dock** on the dialog box and replace the buffer container with one containing the buffer for the experiment. Allow the system to **Prime** by hitting **Start**.
5. **Prime** the system twice by selecting the option under the **Working Tools** in the **Tools Menu**. Repeat this procedure after the first **Prime** has been completed.

Note: Always prime 3x when switching running buffers as this will remove any of the previous buffer from the system and improve results.

### VI. Optimization of pH for Immobilization.

The following procedure is designed to optimize the pH of the immobilization buffer to maximize immobilization to the surface. Efficient immobilization for amine coupling results from the electrostatic pre-concentration of ligand on the surface. To accomplish this pre-concentration, the ligand must be positively charged and thus must be in an immobilization buffer 0.5 – 1 unit in pH below its pI.

A more detailed description of the amine coupling process can be found in the *BIAApplication Handbook*, section 4.2.


1. Dilute the ligand from 1mg/mL to 30ug/mL by the following:
  - a. Place 97uL of each coupling buffer (pH 7.0, pH 6.0, pH 4.8, Coupling buffers 1,2 and 3, respectively) into separate Biacore plastic tubes.
  - b. Add 3uL of ligand to each tube
  - c. Flick the bottom of the tube to mix
2. Place the plastic vials in the rack positions R2D1, R2D2 and R2D3, respectively
3. Place a vial containing Regeneration Solution 3 (10mM glycine, pH 2.0) into rack position R1D3.
4. From the **Run Menu**, select **Run Sensogram...** to begin taking data
  - a. Select Fc2 for both flow and detection and click OK
  - b. Enter a flow rate of 10uL/min
5. Queue up commands in the queue to inject all three ligand solutions by the following:
  - a. Form the **Command Menu**, select **QUICKINJECT** from the menu
    - i. Enter position R2D1
    - ii. Enter 20uL for the injection volume (2min)
    - iii. Click on the start button
  - b. Repeat the same procedure for positions R2D2 and R2D3. The commands will be saved on the queue and performed in sequence
  - c. Queue another **QUICKINJECT** of 20uL of the Regeneration Solution in position R1D3.

Note: The surface has not been activated using NHS/EDC yet so there should be no amine reactive moieties on the surface and thus there should be no coupling of the ligand to the surface.
6. Wait for the injections to finish.

## Biacore Tutorial

4

7. Compare the results with the following sensogram:



**Figure 3.** Anticipated Results from an Immobilization Buffer Optimization

There are several things to note about the sensogram obtained:

- i. At pH 7.0 (the first injection) the pH is very close to the pI of the protein. In these conditions, the protein carries a minimal positive charge and thus a minimal attraction to the negatively charged matrix of the CM5 chip.
  - ii. At pH 6.0 and 4.8, the protein becomes more progressively positively charged which increases the affinity of the analyte for the matrix.
  - iii. The pH 6.0 buffer is the optimum choice as it creates enough positive charge on the protein to attract it to the matrix but at a pH that provides a milder environment for the protein. As a rule of thumb, *the level of attraction should always exceed the desired immobilization level.*
8. Select the **Run Menu** and select **Stop Sensogram** from it. Click **Yes** to stop data collection and **Yes** if the program asks to clear the command queue.
  9. Use the **Save As** function under the File Menu to save the sensogram under a folder with the user's name on it.
  10. Place the instrument in **Standby** by selecting the **Tools Menu**, then **Working Tools** and finally, **Standby**. Ensure enough buffer solution is in the buffer container to last for a few days (~80mL).

### VII. Immobilization and Surface Binding Capacity Check.


1. Mix a fresh 1mL solution of 0.5M NHS/0.2 M EDC. Place 200uL of this solution into a plastic Biacore vial in rack position R2A1.
2. Mix 4uL of the ligand ( $\alpha$ -Myo MAB) stock solution with 396uL of 5mM maleate (pH 6.0, Coupling Buffer 2). Place 100uL of this solution into a Biacore vial in rack position R2A2.
3. Place 200uL of 1M Ethanolamine Hydrochloride (pH 8.5) into a vial in rack position R2A3.
4. Place ~2mL of degassed HBS-EP Buffer into a glass vial in rack position R2F3.

## Biacore Tutorial

5

5. Select the **Run Menu** and select **Start Sensogram...**
  - a. Set the system to flow over Fc1 and detect over Fc1-4
  - b. Set a flow rate of 5uL/min
6. Queue the following commands (from the **Command Menu**):
  - a. **QUICKINJECT** of 35uL of R2A1 (the NHS/EDC solution)
  - b. **QUICKINJECT** of 5uL of R2F2 (HBS-EP)
  - c. **QUICKINJECT** of 35uL of R2A2 (the ligand)
  - d. **QUICKINJECT** of 5uL of R2F2 (HBS-EP)
  - e. **QUICKINJECT** of 35uL of R2A3 (Ethanolamine)
  - f. **QUICKINJECT** of 5uL of R2F2 (HBS-EP)
  - g. **FLOWPATH Fc2** (switches flow to Fc2)
  - h. **QUICKINJECT** of 35uL of R2A1 (the NHS/EDC solution)
  - i. **QUICKINJECT** of 5uL of R2F2 (HBS-EP)
  - j. **QUICKINJECT** of 35uL of R2A3 (Ethanolamine)
  - k. **QUICKINJECT** of 5uL of R2F2 (HBS-EP)

Note: This conditions two parallel surfaces, one containing the ligand for kinetic analysis and another surface with no ligand as a reference.
7. Compare the results gathered with the following figure:



- a. Baseline
- b. EDC/NHS activation
- c. After activation, the baseline is increased by 50-500 RU.
- d. Attraction of ligand into the matrix and covalent coupling.
- e. Buffer washes away loosely associated ligand. The response level at this point (bound) gives a first indication of the immobilized amount.
- f. Deactivation and further washing away of loosely associated ligand.
- g. The difference in response between points 1 and 7 reflects the amount immobilized (final).

**Figure 4.** Sample Immobilization Sensogram.

8. Calculate the theoretical binding capacity of the surface ( $R_{max}$ ) using the following formula:

$$R_{max} = \frac{MW_{analyte}}{MW_{ligand}} * \text{Amount Immobilized} * \text{Stoichiometric Ratio}$$

$MW_{analyte}$  = Molecular Weight of the analyte = 17500 Daltons (myoglobin)

$MW_{ligand}$  = Molecular Weight of the ligand = 156000 Daltons ( $\alpha$ -myoglobin MAB)

Amount Immobilized = RU difference between points a and g

## Biacore Tutorial

Stoichiometric Ratio = 2 (Each antibody can bind two antigens)

The  $R_{\max}$  Obtained from the Sample Sensogram is 1571 RU. Generally, 50-400 RU surface binding capacity is considered optimal for most kinetics experiments.

9. Select the **Run Menu** and select **Stop Sensogram**
10. Save the sensogram as was done before (**File Menu** then **Save As**)
11. Ensure the system is put into **Standby** so the system isn't run dry

### VIII. Regeneration and Nonspecific Binding Test

The purpose of this exercise is to determine the optimum regeneration conditions for the analyte-ligand pair using a series of low pH solutions. In order to retain ligand functionality on the surface, the solutions will be tested from mild to harsh conditions.


An additional step in this analysis is to estimate the non-specific binding of the analyte by the surface. This is an important validation step that verifies that the observed binding phenomena are ligand-analyte interactions and not non-specific adsorption of the analyte.

1. Dilute the analyte solution by adding 24uL of the stock myoglobin solution to 376uL of the HBS-EP buffer. Final concentration should be ~5ug/mL. Cap the vial, mix by flicking the base of the vial and place the vial in R2A1.
2. Place 200uL of Regeneration solutions 1 (10mM acetate, pH 4.0), 2 and 3 (10mM glycine pH 3.0 and 2.0, respectively) into separate Biacore plastic vials and place the vials in rack positions R2B1, R2B2 and R2B3, respectively.
3. Select the **Run Menu** and select **Run Sensogram...**
  - a. Set the detection and flowpath to Fc1-2
  - b. Set the flow rate to 10uL/min
4. Inject the analyte solution by selecting the **Command Menu** and then selecting **Inject...**
  - a. **KINJECT** 30uL of the solution R1A1. Enter a disassociation time of 180s
5. Queue up the following commands:
  - a. **QUICKINJECT** 10uL of R2B1 (Regeneration Solution 1, ph 4.0)
  - b. **QUICKINJECT** 10uL of R2B2 (Regeneration Solution 2, ph 3.0)
  - c. **KINJECT** 30uL of the solution R1A1. The purpose of injecting the analyte again is to ensure that the regeneration doesn't harm the ligand on the surface. Disassociation time again should be 180s
  - d. **QUICKINJECT** 10uL of R2B3 (Regeneration Solution 1, ph 2.0)

## Biacore Tutorial

7

6. Compare the sensogram obtained with the figure below:



### Regeneration Sensogram Evaluation.

- a. The first injection at pH 4.0 shows no regeneration. Bound material remains on the surface after the regeneration pulse. Some dissociation can be seen in the buffer flow after the analyte injection.
- b. The second injection at pH 3.0 shows insufficient regeneration. Part of the bound material remains on the surface after the regeneration pulse.
- c. The third injection at pH 2.0 is good regeneration. All bound material is removed and the baseline is back to about its original level.

**Figure 5:** Sample Regeneration Experiment

7. Next examine the sensogram for Fc2 for non specific binding.
  - a. The baseline RU value before and after injection with the analyte should be equal if no non-specific binding is occurring.
8. The stability of the surface after repeated regeneration must be tested to ensure no ligand is lost with repeated regenerations. Queue up the following commands:
  - a. **KINJECT** 30uL of R1A1 (ligand), Dissociation time: 180s
  - b. **QUICKINJECT** 10uL of R2B3 (Regeneration Solution 1, ph 2.0)
  - c. **KINJECT** 30uL of R1A1 (ligand), Dissociation time: 180s
  - d. **QUICKINJECT** 10uL of R2B3 (Regeneration Solution 1, ph 2.0)
  - e. **KINJECT** 30uL of R1A1 (ligand), Dissociation time: 180s
  - f. **QUICKINJECT** 10uL of R2B3 (Regeneration Solution 1, ph 2.0)
9. Examine the resulting sensogram taking note of the baseline after each regeneration injection. If the baseline decreases after subsequent regeneration cycles, than the regeneration is probably damaging/removing the ligand on the surface and is too harsh. If the baseline slowly rises after

## Biacore Tutorial

---

---

subsequent regeneration cycles, than analyte is being left on the surface and the regeneration needs to be tweaked.

10. Save the sensogram in the proper directory and ensure the system is put into **Standby**.