

The Cardinal Principles of Panel Design


With the arrival of flow cytometers designed to measure 10 or more fluorescent markers, flow cytometry panels require very careful planning and design to avoid some of the inherent problems of spectral overlap. These guidelines provide specific assistance for creating a new or modifying an existing multicolor panel on each of the Cancer Center FCSR flow cytometers. Page 5 also describes basic “backbone panels” of fluorochromes that generally work together well on each instrument.

1. Choose fluorochromes that are spread across the spectrum of lasers and detectors. Selecting the first color on each laser path will create a panel needing very little, if any, compensation. For example, the Yeti has UV, 405, 488, 561, and 604 lasers, so choosing BUV395, BV421, FITC, PE, and APC for a five-color panel is the best choice. This panel will need a very small amount of compensation for PE spillover into FITC (<2%) yielding very little digital spread.
2. Understand the biology of your experiment and define your analysis objectives. Optimized Multicolor Immunofluorescence Panels (OMIP) can be a good place to start. These peer-reviewed panels are available in Cytometry Part A through the Wiley Online Library. The OMIPs provide a wide variety of well-designed analyses for various types of cells, but be sure to verify that color choices are best for your instrument of choice.

Objective	Briefly describe the purpose of your experiment, e.g. “Measuring mouse PBMC subsets”
Species	Human, mouse, etc.
Cell types	List the cell type you will be using, for example PBMC, bone marrow, WBCs in tissue digest

- Sketch a gating scheme to identify marker relationships such as mutually exclusive marker sets and the relationships of low and high expression markers.

Gating Example- OMIP-033: <https://doi.org/10.1002/cyto.a.22889>



- Evaluate the expression of each marker based these criteria (Mahnke YD, *Clin Lab Med.* 2007; 27:469-485)

Primary- well characterized expression which is positive or negative.


Secondary- mid-range expression or a continuum of expression.

Tertiary- expressed at low levels, unknown levels, or is rare cell type.

Manufacturers frequently provide expression levels on their websites. BioLegend is particularly helpful. www.biolegend.com

- List the color choices available for each marker. Include what your lab has in stock, but be aware that these may not be the best choices for a large parameter panel. Making a table is useful here.

- Evaluate the stain index of the available fluorochromes. Note that stain indices can vary between instrument models. High (primary) expression markers should be used with low stain index flours and tertiary (low) expression should be on flours with a high stain index.




- Evaluate the spillover behavior of each available fluorochrome using the Spillover Spreading Matrix (SSM). The SSM summarizes how flours behave together by measuring the amount of signal spread into the other detector channels. The flours in the Y column spread into the channels of the X columns. While SSM is directly related to the amount of compensation, it is independent of signal intensity. **Fluorochromes pairs that give and receive high spillover should be used for mutually exclusive markers. Low (tertiary) expression markers should be on flours that don't receive spread.**
- Use this information to create a table and assign color choices. See Table I.

Table I: Example of summary information and fluorochrome choices for a simple panel.

Laser, Filter, & Common flours	Marker	Flour	Stain Index	Purpose	Expression Primary Secondary Tertiary
387/11 BUV395	CD3	BUV395	high	T cells	primary
420/10 BV421	CD19	BV421	high	B cells	primary
509/24 FITC/GFP	CD45	FITC	low	WBC	primary
577/15 PE	CD34	PE	high	Hematopoetic progenitors	Tertiary (rare)
670/30 APC/A647	EpCAM	Alexa 647	high	epithelial	primary
775/50 APC Cy7	Live-dead	NIR	low	Live-dead	Secondary

- At this point, we would be pleased to evaluate your panel. Please bring all of the above information, including the purpose, gating scheme, available reagents, and the summary information table.
- Test your panel. Many manufacturers sell trial samples of antibodies for testing. Be sure to titer all antibodies and include FMO controls to delineate negative from positive. **Choosing the titer that gives a good separation rather than the highest signal/noise ratio will reduce spreading into other fluorochromes**



Panel Design Cuidelines
 Created by CC Flow Cytometry Shared Resource
 Date Created: 9/26/2019
 Date Modified: 9/26/2019

Backbone panels of fluorochromes that generally work together well:

Gallios Backbone Panels

Laser (Filter)	4 Color Panel	8 Color Panel	Notes
488 FL1 (525/40)	FITC or Alexa 488	FITC or Alexa 488	
488 FL2 (575/30)	PE	PE	
488 FL3 (620/30)			
488 FL4 (695/30)		PE Cy5.5	
488 FL5 (755LP)		PE Cy7	
638 FL6 (660/20)	APC or Alexa 647	APC or Alexa 647	
638 FL7 (725/20)			
638 FL8 (755LP)		APC Cy7	Live/dead NIR, Zombie NIR
405 FL9 (450/40)	BV421 or Pacific Blue	BV421 or Pacific Blue	
405 FL10 (550/40)		BV510	

Yeti Backbone Panels

Laser (Filter)	6 Color Panel	10 Color Panel	Notes
355(387/11)	BUV395	BUV395	
405(420/10)	BV421	BV421	
488(509/24)	FITC	FITC	
561(577/15)	PE	PE	
641(670/30)	APC	APC	
641(775/50)	APC Cy7	APC Cy7	Live/dead NIR, Zombie NIR
405(525/59)		BV510	
405(750LP)		BV 711	
561(750LP)		PE Cy7	
488(692/80)		PerCP cy5.5	

Aurora Backbone 24 Color Panel (careful planning required)


Violet		Blue	Yellow	Red
BV421	BV650	BB515	PE	APC
Super Bright 436	BV711	AF488	PE DAZZLE 594	AF647
eFluor450	BV750	AF532	PECy5	A700
BV480	BV785	PerCP Cy5.5	PECy7	APC Fire 750
BV510		PerCP eFluor 710		
BV570				
BV605				

Panel Design Cuidelines
 Created by CC Flow Cytometry Shared Resource
 Date Created: 9/26/2019
 Date Modified: 9/26/2019

Yeti Experiment Summary Table

Laser, Filters, & fluors	Marker	Fluor	Stain Index	Purpose	Expression
387/11 BUV395					
447/60 BUV496					
525/50 BUV563					
670/30 BUV661					
700LP BUV737					
420/10 BV421					
460/22 PACIFIC BLUE					
525/50 BV510					
615/24 BV650					
670/30 BV650					
720/60 BV711					
750LP BV786					
509/24 FITC/GFP					
549/15 YFP					
692/80 PerCP					
750LP PE Cy7**					
577/15 PE					
589/15 tdTomato					
615/24 PE Dazzle					
640/20 mPlum					
670/30 PE Cy5					
720/60 PE Cy5.5					
750LP PE Cy7					
670/30 APC/A647					
720/60 A700					
775/50 APC Cy7					
800LP APC ef780					

Yeti Stain Index




Yeti Spillover Spreading Matrix

	FL3	FL5	FL6	FL7	FL8	FL11	FL13	FL14	FL15	FL16	FL17	FL18	FL19	FL20	FL26	FL27	FL29
	FITC	PECy7(488)	PCP-Cy5.5	PECY7	PECY5	PE	PE-Dazzle	BV650	BV711	BV786	PacBl	BV421	BV605	BV510	Alexa700	APC Fire75	Alexa 647
FITC	1.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
PECy7	0.0	1.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.2	0.0	0.0	0.0	0.0	0.0	0.1	0.0
PCP-Cy5.5	0.0	0.9	1.0	0.4	0.4	0.0	0.0	0.2	1.2	0.7	0.0	0.0	0.0	0.0	0.8	0.2	0.2
PECY5	0.0	0.2	0.7	0.3	1.0	0.0	0.0	0.1	0.9	0.0	0.0	0.0	0.0	0.0	0.3	0.1	0.4
PE	0.0	0.0	0.1	0.0	0.1	1.0	0.8	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
PE-Dazzle	0.0	0.0	0.1	0.0	0.2	0.2	1.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
BV650	0.0	0.0	0.0	0.0	0.1	0.0	0.1	1.0	1.4	0.6	0.0	0.1	0.1	0.0	0.2	0.0	0.2
BV711	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	1.0	0.6	0.0	0.0	0.0	0.0	0.2	0.1	0.0
BV786	0.0	0.1	0.0	0.0	0.0	0.0	0.0	0.0	0.3	1.0	0.00322	0.0	0.0	0.0	0.0	0.0	0.0
PacBl	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	1604.0	0.0	1.0	0.1	0.0	0.5	0.0	0.0	0.0
BV421	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.3	1.0	0.0	0.1	0.0	0.0	0.0
BV605	0.0	0.0	0.1	0.0	0.3	1.7	2.3	0.8	1.2	0.5	0.3	0.8	1.0	0.1	0.0	0.0	0.0
BV510	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.1	0.2	0.1	0.1	0.0	0.1	1.0	0.0	0.0	0.0
Alexa700	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.1	0.1	0.0	0.0	0.0	0.0	1.0	0.3	0.0
APC Fire750	0.0	0.0	0.0	0.2	0.0	0.0	0.0	0.0	0.0	0.2	0.0	0.0	0.0	0.0	0.0	1.0	0.0
Alexa 647	0.0	0.0	0.0	0.0	0.2	0.0	0.0	4812.0	0.0	0.0	0.0	0.0	0.0	0.0	0.7	0.2	1.0

Panel Design Cuidelines
 Created by CC Flow Cytometry Shared Resource
 Date Created: 9/26/2019
 Date Modified: 9/26/2019

Gallios Experiment Design Summary

Laser, Filters& Common fluors	Marker	Fluor	Stain Index	Purpose	Expressio n
488 FL1 (525/40)	FITC or Alexa 488				
488 FL2 (575/30)	PE				
488 FL3 (620/30)	PE-Texas Red, ECD				
488 FL4 (695/30)	PE Cy5 or PE Cy5.5				
488 FL5 (755LP)	PE Cy7				
638 FL6 (660/20)	APC or Alexa 647				
638 FL7 (725/20)	Alexa 700				
638 FL8 (755LP)	APC Cy7				
405 FL9 (450/40)	BV421 or Pacific Blue				
405 FL10 (550/40)	BV510				



Gallios Spillover Spreading Matrix

SSM	FL1 FITC	FL2 PE	FL3 PE Dazzle	FL4 PE Cy5.5	FL5 PE Cy7	FL6 APC	FL7 A700	FL8 APC FIRE750	FL9 BV421	FL10 BV510
FL1 FITC	1	.125.4	0.08938	0.01229	0.00288	0	0	0.0003	0.00106	0.07678
FL2 PE	0.0123	1	0.83608	0.13055	0.02174	0	0	0	0.00059	0.19051
FL3 PE Dazzle	0.00231	0.17172	1	0.2096	0.03975	0.00101	0.00036	0	0.00146	0.03187
FL4 PE Cy5.5	0	0	0	1	0.33902	0.0515	0.32253	0.10127	0.00193	0
FL5 PE Cy7	0.00039	0.00343	0.00291	0.00118	1	0	0.0078	0.11654	0.00091	0.00084
FL6 APC	0	0	0	0.00227	0.00086	1	0.32144	0.07461	0.00129	0.00159
FL7 A700	0	0	0	0.00443	0.00493	0.00288	1	0.22571	0.00393	0.00272
FL8 APC FIRE750	0	0	0	0	0.01651	0.03764	0.06476	1	0.00552	0.00364
FL9 BV421	0	0	0	0	0.00031	0	0	0	1	0.11775
FL10 BV510	0.0001	0	0	0.0001	0.00031	0	0	0.00005	0.07195	1

Panel Design Cuidelines
Created by CC Flow Cytometry Shared Resource
Date Created: 9/26/2019
Date Modified: 9/26/2019

Cytek Aurora Panel Design

The spectral fluorochrome signatures, lists of unique versus similar dyes, stain index tables, and spillover spreading matrices for the Aurora are available on the Cytek Biosciences website in the Data Sheets section of the Resources tab.

<https://cytekbio.com/blogs/resources/tagged/data-sheets>

All of the previously discussed Cardinal Rules of Panel Design apply to panels for the Aurora and spreading can still be a major problem. Careful attention to the details is very important. Panel evaluation should include inspecting the unmixing of each reference control by analyzing each control against each of the other colors in the panel. Build a worksheet with a dual-parameter plot for each of the fluorochromes, placing each fluor on a y-axis. Inspect each reference control by assigning that control fluorochrome to all of the x axes. Look for swoops and other odd angles that indicate over or under correction of the spillover. Slight corrections can be made by using the manual compensation option in the Aurora software.

For very large Aurora panels, breakdown the panel into plots of mini-panels and inspect the plots to verify that the data is biologically correct. For example, look at all of the T cell plots and examine them for odd populations which do not make biological sense or have other inconsistencies.