

# Duolink<sup>®</sup> In Situ

## SHORT INSTRUCTIONS – FLUORESCENCE

### 1. BLOCKING

Add blocking solution to samples.  
Incubate.



▶ *Remove block*

### 2. PRIMARY ANTIBODIES

Dilute the primary antibodies in appropriate buffer and apply to samples.  
Incubate.



▶ *Wash in suitable buffer for 2 × 5 min*

### 3. PLA<sup>®</sup> PROBES

Dilute the two PLA probes 1:5 in appropriate buffer and apply to samples.  
Incubate for 60 min at +37°C.



▶ *Wash in 1x Wash Buffer A for 2 × 5 min*

### 4. LIGATION


Dilute the Ligation stock 1:5 in H<sub>2</sub>O. Dilute the Ligase at 1:40 in the solution and apply the mix to samples.  
Incubate for 30 min at +37°C.



▶ *Wash in 1x Wash Buffer A for 2 × 2 min*

### 5. AMPLIFICATION

Dilute the Amplification stock 1:5 in H<sub>2</sub>O. Dilute the Polymerase at 1:80 in the solution and apply the mix to samples.  
Incubate for 100 min at +37°C.




▶ *Wash in 1x Wash Buffer B for 2 × 10 min*

▶ *Wash in 0.01x Wash Buffer B for 1 min*

### 6. PREPARATION FOR IMAGING

Mount the samples with Duolink In Situ Mounting Medium with DAPI,  
wait for 15 min, and analyze in a fluorescence or confocal microscope.



*Note: Use open droplet reactions without a cover slip and perform all incubations in a humidity chamber. Use volumes corresponding to your reaction area, see the Reaction Volume Guide (# 0062). Use a freezing block when removing the enzymes from the freezer (-20°C). Washing should be done in a minimum volume of 70 ml on a shaker with gentle orbital shaking.*