

Attune NxT Parameters List

Lasers	Filter	Parameter	Common Fluorophores
V405	450/40	VL1	BV421 /Pacific Blue /E450 /A405 /Zombie Violet /V450 /BFP, mCerulean /Azurite /GhostViolet450 /DAPI
	525/50	VL2	BV510 /Zombie Aqua /V500 /A430 /E506 /Sapphire /AmCyan /GhostViolet510 /GhostViolet540
	610/20	VL3	BV605
	660/20	VL4	BV650
	710/50	VL5	BV711
	780/60	VL6	BV786
B488	FSC 488/10	BL FSC	
	SS 488/10	BL SS	
	530/30	BL1	FITC /A488 /GFP /YFP /Zombie Green /BB515 /Dylight488 /ZSGreen /Emerald /AzamiGreen /mWasabi /Midori-Ishi /Cyan /GhostBlue516 /A532 /E520 /Dylight549 /Rhodamine123 /Venus /Topaz /mCitrine
	695/40	BL2	BB700 /PerCPcy5.5 /PerCpe710
Y561	586/16	YLI	PE /dTomato /A546 /A555 /mOrange2 mBanana /KusabiraOrange /mApple /A568 /E570 /Dylight594 /RFP /DsRed /mTangerine /mStrawberry
	620/15	YL2	PE-C594 /A594 /PEe610 /e615 /mCherry /mRuby /PE-Dazzle /mRFP/mRaspberry /ZombieRed /PI
	780/60	YL3	PE-Cy7
R637	670/14	RLI	APC /A647 /A660 /E660 /Dylight649 /RFP657
	720/30	RL2	A700 /APC-700 /mIFP /ZombieNIR /GhostRed710
	780/60	RL3	Live/Dead NIR /APCy7 /APCh7 /APCe780 /A790 /GhostRed780